

STANDARD OPERATING PROCEDURE

SOP # 3 Page: 1 of 2	Version 1.1 Date: 10/08/2017 Revise by: 10/08/2017	Prepared By: Dr Joyanta Modak Approved By: Dr Vijay Dhanasekaran
-----------------------------	--	---

Laboratory of Viral Evolution, Department of Microbiology, Monash University

SOP 3: Fecal sample process

1. INTRODUCTION

This SOP describes the procedure for fecal sample processing from different animals (e.g. Bats, Possums, dogs, cats, ducks, swamphens).

2. SCOPE

This procedure serves as a guide to laboratory personnel engaged in fecal sample processing from different sources.

3. SAFETY

All lab workers must:

- Wear appropriate personal protective equipment (PPE) (Nitrile gloves, lab coat, closed-toed shoes, and safety glasses when specified by risk assessment) when carrying out the procedure.
- Sample handling has to be carried out in PC2 facilities.
- USE RNaseZap spray to remove RNase contamination from work surface/ equipment/ solutions.
- Use disposable, individually wrapped, sterile plastic ware and sterile, disposable RNase-free pipettes, pipette tips, and tubes.
- Read and understand the recommendations in this SOP to ensure uniformity in practice.
- New staff or students should seek supervision from a senior lab co-worker during first attempt of the procedure.

4. Materials Needed

- Hanks' balanced salt solution (Gibco BRL)
- 0.45-µm filter and a 0.22-µm filter (Millipore)
- Centrifuge that reach to 12000 g at 4°C
- Ultracentrifuge 184,000 × g
- Turbo DNase (final concentration, 20 U/ml; Ambion)
- RNase A (final concentration, 0.1 mg/ml; Fermentas)
- Falcon tube

STANDARD OPERATING PROCEDURE

SOP # 3 Page: 2 of 2	Version 1.1 Date: 10/08/2017 Revise by: 10/08/2017	Prepared By: Dr Joyanta Modak Approved By: Dr Vijay Dhanasekaran
-----------------------------	--	---

Laboratory of Viral Evolution, Department of Microbiology, Monash University

5. Fecal Sample processing

The fecal samples will be transported to the laboratory in dry ice or liquid nitrogen and stored at -80°C until analysis.

1. Resuspend the fecal sample in Hanks' balanced salt solution (Gibco BRL) (1:10, wt/vol) and homogenized thoroughly.
2. Centrifuge the suspension at $5,000 \times g$ for 10 min.
3. Transfer the supernatant to fresh tubes and centrifuged at $13,000 \times g$ for 15 min.
4. Filter a volume of 150 ml of the supernatant through a $0.45\text{-}\mu\text{m}$ filter and a $0.22\text{-}\mu\text{m}$ filter (Millipore) to remove eukaryotic and bacterial cell-sized particles.
5. Ultracentrifuge the filtrate at $184,000 \times g$ for 3 h at 8°C to pellet the viral particles.
6. Resuspend the pellets in $400\text{ }\mu\text{l}$ Hanks' solution and filter with a $0.22\text{-}\mu\text{m}$ filter.
7. Treat the filtrate with Turbo DNase (final concentration, 20 U/ml; Ambion) and RNase A (final concentration, 0.1 mg/ml; Fermentas) at 37°C for 30 min to digest non-particle-protected nucleic acids.
8. Use the treated sample for RNA extraction.